

Experiment no : 02

Experiment : Determination of the extinction co-efficient of Cupra – ammonium complex , $(\text{Cu}(\text{NH}_3)_2^{2+})$ and determination of the concentration of given unknown solution.

Theory:

According to Beer-Lambert's law.

$$A = \epsilon Cl$$

Where;

A=Absorbance (Optical density)

C= concentration of the solution

L= Light path length

ϵ =Extinction coefficient

Procedure:

- 1) Firstly, prepare a 0.1 mol dm^{-3} CuSO_4 solution in a 100 mL volumetric flask. Transfer about 10.00 mL of 0.1 mol dm^{-3} CuSO_4 solution into 50 mL volumetric flask and add 10.00 mL of conc NH_3 . (Resulting solution is called as Cupra-ammonium stock solution).
- 2) Take 10.00 mL, 8.00 mL, 6.00 mL, 5.00 mL, 4.00 mL and 2.00 mL of the above stock solution into 50.00 mL volumetric flasks using a burette and dilute them to 50.00 mL with distilled water.
- 3) Determine the best filter for the Cupra-ammonium complex by measuring the absorbances of lowest and highest concentrations at each wave length.
- 4) Measure the absorbances of the above standard series solutions at the appropriate best filter.
- 5) Plot absorbance vs concentration (mol dm^{-3}) of Cu – ammonium complex .
- 6) Determine the extinction coefficient of the Cupra –ammonium complex. Find the concentration of an unknown Cu^{2+} solution. (Take only 10.00 mL of the unknown Cu-ammonium complex solution)