


Lecture 02



INFRARED SPECTROSCOPY



THE ELECTROMAGNETIC SPECTRUM




- Infrared radiation stimulates molecular vibrations.
- Infrared spectra are traditionally displayed as %T (percent transmittance) versus wave number (4000-400 cm^{-1}).
- Useful in identifying presence or absence of functional groups.



Wavenumber = $1/\text{wavelength}$


INFRARED SPECTROSCOPY

- In the IR region of the electromagnetic spectrum, the absorption of radiation by a sample is due to changes in the vibrational energy states of a molecule.




INFRARED SPECTROSCOPY


Methane




Asymmetrical stretching




Symmetrical stretching




Bending or scissoring



Rocking or in plane bending



Twisting or out-of-plane bending




Wagging or out-of-plane bending


General Theory of IR Spectroscopy

Normal modes of vibration = stretching and bending

Stretching vibrations




symmetric stretch



asymmetric stretch


Bending vibrations




symmetric in-plane
bend (scissor)



asymmetric in-plane
bend (rock)




symmetric out-of-plane
bend (twist)



asymmetric out-of-plane
bend (wag)

INFRARED SPECTROSCOPY

- Only vibrations that cause a change in 'polarity' give rise to bands in IR spectra



Symmetric stretch

Asymmetric stretch


Bending (doubly degenerate)

INFRARED SPECTROSCOPY

- Any change in shape of the molecule- stretching of bonds, bending of bonds, or internal rotation around single bonds
- Asymmetrical stretching/bending and internal rotation change the dipole moment of a molecule. Asymmetrical stretching/bending are IR active.
- Symmetrical stretching/bending does not. Not IR active

INFRARED SPECTROSCOPY


- Human Breath



INFRARED SPECTROSCOPY


- How much movement occurs in the vibration of a C-C bond?

stretching vibration



For a C-C bond with a bond length of 154 pm, the variation is about 10 pm.



bending vibration




For C-C-C bond angle a change of 4° is typical. This moves a carbon atom about 10 pm.

INFRARED SPECTROSCOPY


- Infrared (IR) electromagnetic radiation causes vibrations in molecules (wavelengths of 2500-15,000 nm or 2.5 – 15 mm)
- The greater the mass - the lower the wavenumber



INFRARED SPECTROSCOPY



Ethane




Chloroethane

INFORMATION OBTAINED FROM IR SPECTRA

- IR is most useful in providing information about the presence or absence of specific **functional groups**.
- IR can provide a **molecular fingerprint** that can be used when comparing samples. If two pure samples display the same IR spectrum it can be argued that they are the same compound.
- IR **does not** provide detailed information or proof of molecular formula or structure. It provides information on molecular fragments, specifically functional groups.
- Therefore it is very limited in scope, and must be used in conjunction with other techniques to provide a more complete picture of the molecular structure.


IR ABSORPTION RANGE

- The typical IR absorption range for covalent bonds is **600 - 4000 cm^{-1}** . The graph shows the regions of the spectrum where the following types of bonds normally absorb. For example a sharp band around 2200-2400 cm^{-1} would indicate the possible presence of a C-N or a C-C triple bond.



THE FINGERPRINT REGION

- Although the entire IR spectrum can be used as a fingerprint for the purposes of comparing molecules, the **600 - 1400 cm^{-1}** range is called the **fingerprint region**. This is normally a complex area showing many bands, frequently overlapping each other. This complexity limits its use to that of a fingerprint, thus neglect the region left of 1400 cm^{-1} .



Focus your analysis on this region. This is where most stretching frequencies appear.


Fingerprint region: complex and difficult to interpret reliably.

Characteristic IR frequencies of common functional groups

Bond	Type of Compound	Frequency Range, cm^{-1}	Intensity
C—H	Alkanes	2850–2970 1340–1470	Strong Strong
C—H	Alkenes (>C=C<H)	3010–3095 675–995	Medium Strong
C—H	Alkynes ($\text{—C}\equiv\text{C—H}$)	3300	Strong
C—H	Aromatic rings	3010–3100 690–900	Medium Strong
O—H	Monomeric alcohols, phenols	3590–3650	Variable
	Hydrogen-bonded alcohols, phenols	3200–3600	Variable, sometimes broad
	Monomeric carboxylic acids	3500–3650	Medium
	Hydrogen-bonded carboxylic acids	2500–2700	Broad
N—H	Amines, amides	3300–3500	Medium
C=C	Alkenes	1610–1680	Variable
C=C	Aromatic rings	1500–1600	Variable
C $\equiv$ C	Alkynes	2100–2260	Variable
C—N	Amines, amides	1180–1360	Strong
C $\equiv$ N	Nitriles	2210–2280	Strong
C—O	Alcohols, ethers, carboxylic acids, esters	1050–1300	Strong
C=O	Aldehydes, ketones, carboxylic acids, esters	1690–1760	Strong
NO ₂	Nitro compounds	1500–1570 1300–1370	Strong Strong


IR SPECTRUM OF ALKANES

- Alkanes have no functional groups. Their IR spectrum displays only C-C and C-H bond vibrations. Of these the most useful are the **C-H bands**, which appear around **3000 cm^{-1}** . Since most organic molecules have such bonds, most organic molecules will display those bands in their spectrum.




IR SPECTRUM OF ALKENES

- Besides the presence of C-H bonds, alkenes also show sharp, medium bands corresponding to the **C=C bond stretching vibration** at about **1600-1700 cm⁻¹**. Some alkenes might also show a band for the =C-H bond stretch, appearing around **3080 cm⁻¹** as shown below. However, this band could be obscured by the broader bands appearing around 3000 cm⁻¹



IR SPECTRUM OF ALKENES


- This spectrum shows that the band appearing around **3080 cm^{-1}** can be obscured by the broader bands appearing around **3000 cm^{-1}** .



IR SPECTRUM OF ALKYNES


- The most prominent band in alkynes corresponds to the **carbon-carbon triple bond**.
- It shows as a sharp, weak band at about **2100 cm^{-1}** . The reason it's weak is because the triple bond is not very polar.
- In some cases, such as in highly symmetrical alkynes, it may not show at all due to the low polarity of the triple bond associated with those alkynes.
- Terminal alkynes have C-H bonds involving the *sp* carbon (the carbon that forms part of the triple bond). Therefore they may also show a sharp, weak band at about **3300 cm^{-1}** corresponding to the C-H stretch.
- Internal alkynes, that is those where the triple bond is in the middle of a carbon chain, do not have C-H bonds to the *sp* carbon and therefore lack the aforementioned band.

IR SPECTRUM OF ALKYNES



IR SPECTRUM OF A NITRILE

In a manner very similar to alkynes, nitriles show a prominent band around **2250 cm^{-1}** caused by the **CN triple bond**. This band has a sharp, pointed shape just like the alkyne C-C triple bond, but because the CN triple bond is more polar, this band is stronger than in alkynes.



Signal Intensity in an IR Spectrum


For a bond to absorb in IR, there must be a change in the dipole moment during the vibration.

If a bond is non-polar due to symmetry, IR light will not be absorbed and the vibration is said to be IR inactive.



$\text{C}\equiv\text{C}$ is IR active


symmetry



$\text{C}\equiv\text{C}$ is IR inactive

IR SPECTRUM OF AN ALCOHOL


- The most prominent band in alcohols is due to the **O-H bond**, and it appears as a strong, broad band covering the range of about **3000 - 3700 cm^{-1}** . The sheer size and broad shape of the band dominate the IR spectrum and make it hard to miss.



IR SPECTRUM OF ALDEHYDES AND KETONES


- Carbonyl compounds **contain C=O functional group**.
- Aldehydes and ketones show a strong, prominent, stake-shaped band around **1710-1720 cm^{-1}** (right in the middle of the spectrum). This band is due to the **highly polar C=O bond**. Because of its position, shape, and size, it is hard to miss.
- Because aldehydes also contain a C-H bond to the sp^2 carbon of the C=O bond, they also show a pair of medium strength bands positioned about **2700 and 2800 cm^{-1}** . These bands are missing in the spectrum of a ketone because the sp^2 carbon of the ketone lacks the C-H bond.

IR SPECTRUM OF ALDEHYDES AND KETONES



IR SPECTRUM OF A CARBOXYLIC ACID


- A carboxylic acid functional group combines the features of alcohols and ketones because it has both the **O-H bond** and the **C=O bond**. Therefore carboxylic acids show a very strong and broad band covering a wide range between **2800** and **3500 cm^{-1}** for the O-H stretch. At the same time they also show the stake-shaped band in the middle of the spectrum around **1710 cm^{-1}** corresponding to the C=O stretch.



IR SPECTRA OF AMINES


The most characteristic band in amines is due to the **N-H bond stretch**, and it appears as a weak to medium, somewhat broad band (but not as broad as the O-H band of alcohols). This band is positioned at the left end of the spectrum, in the range of about **3200 - 3600 cm^{-1}** .


Primary amines have two N-H bonds, therefore they typically show two spikes that make this band resemble a molar tooth. Secondary amines have only one N-H bond, which makes them show only one spike, resembling a canine tooth. Finally, tertiary amines have no N-H bonds, and therefore this band is absent from the IR spectrum altogether. The spectrum below shows a secondary amine.



IR SPECTRUM OF AMIDES

The amide functional group combines the features of amines and ketones because it has both the **N-H bond** and the **C=O bond**. Therefore amides show a very strong, somewhat broad band at the left end of the spectrum, in the range between **3100** and **3500 cm^{-1}** for the N-H stretch. At the same time they also show the stake-shaped band in the middle of the spectrum around **1710 cm^{-1}** for the C=O stretch. As with amines, primary amides show two spikes, whereas secondary amides show only one spike.





Benzoic acid

INFRARED SPECTROSCOPY

POSITION	REDUCED MASS BOND STRENGTH (STIFFNESS)	LIGHT ATOMS HIGH FREQUENCY STRONG BONDS HIGH FREQUENCY
STRENGTH	CHANGE IN 'POLARITY'	STRONGLY POLAR BONDS GIVE INTENSE BANDS
WIDTH	HYDROGEN BONDING	STRONG HYDROGEN BONDING GIVES BROAD BANDS

INFRARED SPECTROSCOPY

- In general

Bond →	C-H	C-D	C-O	C-Cl
ν/cm^{-1} →	3000	2200	1100	700
Bond →	C≡O	C=O	C-O	
ν/cm^{-1} →	2143	1715	1100	

INFRARED SPECTROSCOPY

4000-3000 cm⁻¹	3000-2000 cm⁻¹	2000-1500 cm⁻¹	1500-1000 cm⁻¹
O-H N-H C-H	C≡C C≡N	C=C C=O	C-O C-F C-Cl deformations

Increasing energy



Increasing frequency



INTERPRETATION OF INFRARED SPECTRA


- An element of judgement is required in interpreting IR spectra but you should find that it becomes relatively straightforward with practice.
- It is often possible to assign the peaks in the 1600-3600 cm^{-1} region by consulting tables or databases of IR spectra. When making an assignment, give both the type of bond and the type of vibration, *e.g.* O-H stretch or C-H bending vibration.
- The most useful regions are as follows:
 - 1680-1750 cm^{-1} : C=O stretches feature very strongly in IR spectra and the type of carbonyl group can be determined from the exact position of the peak.
 - 2700-3100 cm^{-1} : different types of C-H stretching vibrations.
 - 3200-3700 cm^{-1} : various types of O-H and N-H stretching vibrations.
- Too many bonds absorb in the region of 600-1600 cm^{-1} to allow confident assignment of individual bands. However, this region is useful as a fingerprint of a molecule, *i.e.* if the spectrum is almost identical to an authentic reference spectrum then the structure can be assigned with some confidence.


Infrared Instrumentation

IR Source

Sample
compartment

Detector






- **Dispersive instruments:** with a monochromator to be used in the mid-IR region for spectral scanning and quantitative analysis.
- **Fourier transform IR (FTIR) systems:** widely applied and quite popular in the far-IR and mid-IR spectrometry.
- **Nondispersive instruments:** use filters for wavelength selection or an infrared-absorbing gas in the detection system for the analysis of gas at specific wavelength.

Infrared Instrumentation

- All modern instruments are Fourier Transform instruments.
- In all transmission experiments radiation from a source is directed through the sample to a detector.
- The measurement of the type and amount of light transmitted by the sample gives information about the structure of the molecules comprising the sample.

Infrared Instrumentation



Infrared Experimental

- To obtain an IR spectrum, the sample must be placed in a “container” or cell that is transparent in the IR region of the spectrum.
- Sodium chloride or salt plates are a common means of placing the sample in the light beam of the instrument.

InfraredExperimental



IR transparent Salt Plates

Infrared Experimental

- These plates are made of salt and must be stored in a water free environment such as the dessicator shown here.



Infrared Experimental

- The plates must also be handled with gloves to avoid contact of the plate with moisture from one's hands.



Infrared Experimental

- To run an IR spectrum of a liquid sample, a drop or two of the liquid sample is applied to a salt plate.
- A second salt plate is placed on top of the first one such that the liquid forms a thin film “sandwiched” between the two plates.



Infrared Experimental

- The cell holder is then placed in the beam of the instrument.
- The sample is then scanned by the instrument utilizing predestinated parameters.
- A satisfactory spectrum has well defined peaks-but not so intense as to cause flattening on the bottom of the peaks.



Infrared Experimental

- The salt plates are cleaned by rinsing into a waste container with a suitable organic solvent-commonly cyclohexane - **NEVER WATER!**
- Cloudy plates must be polished to return them to a transparent condition.
- To polish cloudy windows, rotate salt plate on polishing cloth.

